

QUALITY

- High-quality raw HA material
- Improved product viscosity
- Physiologically balanced pH and osmotic pressure for optimal biocompatibility and patient comfort

SAFETY

- Elimination of residual BDDE
- Class-leading glass syringe minimizes the risk of plastic contamination in the HA
- Near-zero endotoxin levels for increased patient safety

EASY TO USE

- Consistent extrusion force from a glass syringe for smooth and controlled injection
- Uniform product delivery to ensure precise and predictable results
- Formulated with lidocaine to enhance patient comfort during treatment



Imported and Distributed by:
PT ADYA ARTHA ABADI

Komplek Pergudangan Multiguna,
Jalan Raya Serpong Blok DI No. 5,
Serpong Utara, Tangerang Selatan,
Banten, Indonesia.

Sales Office:
BERNO AESTHETICS

Ruko Puri Mansion, Blok C-15,
Kembangan, Jakarta Barat,
Indonesia.

bernoaesthetics.com 
info@bernoaesthetics.com 
bernoaesthetics 
+62813 9555 938 



Berno Aesthetics

UTH



INTRODUCING A GROUNDBREAKING RANGE
OF MONOPHASIC DERMAL FILLERS FROM UTH



Why UTH HA Filler?

UTH HA filler is a **cross-linked hyaluronic acid dermal filler** formulated for facial rejuvenation, contouring, and volume restoration.

Each syringe contains **24 mg/ml of high-quality HA** to deliver long-lasting volume and hydration, combined with **0.3% lidocaine** to minimize discomfort during treatment.



	Hydrate 01	Fine 02	Deep 03	Sub-Q 04
Volume	1 ml	1 ml	1 ml	1 ml
Needle	2 x 30G	2 x 30G	2 x 27G	2 x 25G

